


Doppler & Suction



Features:

- Large Colour TFT Screen
- Colour Display with Both FHR & Curve Display
- High-fidelity, crystal clear sound
- High sensitivity Doppler probe
- 2.5Mhz High frequency probe
- Accurate FHR with Loud Volume
- Re-Chargable Battery with Samsung Cord for easy charge
- Auto Shut Off
- Applications: At Home, clinics, community and hospital

Technical Specification:

- FHR Range: 50*240bpm
- Resolution:*1bpm



PORTABLE SUCTION



DOUBLE JAR SUCTION



We support service...